

1 Optical performance of a monofocal intraocular lens designed to
2 extend depth of focus.

3 Fidel Vega, PhD [1]; María S. Millán, PhD [1]; Miguel A. Gil, MD [2]; Nuria
4 Garzón, OD, PhD [3]

5 [1]: Departament d'Òptica i Optometria, Universitat Politècnica de Catalunya,
6 BARCELONATECH, Spain

7 [2]: Hospital de la Santa Creu i Sant Pau, Universitat Autònoma de Barcelona,
8 Spain

9 [3]: Miranza IOA Madrid, Spain

10

11 Corresponding author:

12 Fidel Vega

13 Departament d'Òptica i Optometria. Violinista Vellsolà 37, 08222-Terrassa,
14 Spain

15 Phone:+34 937398333

16 Fax: +34 937398301

17 E-Mail: fidel.vega@upc.edu

18

19 None of the authors have any conflict of interest on the products mentioned in
20 this work

21

22 Supported by project DPI2016-76019-R from the Spanish Ministerio de
23 Economía y Competitividad y Fondos FEDER.

24 **ABSTRACT**

25 **Purpose**

26 To test the performance of a new monofocal intraocular lens, intended to extend
27 depth of focus (Tecnis® Eyhance, ICB00) (ICB-IOL), in comparison to a time-
28 tested standard monofocal IOL (Tecnis® 1-piece, ZCB00) (ZCB-IOL) of same
29 platform and material.

30 **Methods**

31 Assessment of the optical performance of the two IOLs was made *in-vitro* using
32 an optical test bench with a model eye. The spherical aberration (SA),
33 modulation transfer function (MTF) and the area under the MTF (MTFa) were
34 obtained for pupil sizes ranging from 2.0mm to 5.0mm. Through-focus MTFa
35 curves between -3.0D to +1.0D were obtained with three pupils (2.0mm, 3.0mm
36 and 4.5 mm). Halo formation was also assessed for both lenses.

37 **Results**

38 The ICB-IOL had slightly worse optical quality at its best focus (i.e., lower MTF
39 scores at distance vision) and more negative SA than the ZCB-IOL for pupils
40 ranging from 2.0mm up to 3.0mm. The maximum of the through-focus MTFa
41 curve of the ICB-IOL with a 2.0mm pupil, shifted to myopic defocus of -0.50 D.
42 For larger pupils (≥ 3.5 mm), there were no differences of SA, MTF scores and
43 halo energy between the two lenses.

44 **Conclusions**

45 The new ICB-IOL is a modified monofocal lens with 0.50D of additional power in
46 its central 2mm zone and more negative SA values, which induces a myopic
47 shift of the maximum of optical quality and could improve intermediate vision.
48 For pupils larger than 3.5mm, there were no differences between both IOLs.

49 The new ICB-IOL design would produce photic phenomena comparable to a
50 standard IOL.

51

52 **INTRODUCTION**

53 Restoring vision by replacing the opacified crystalline lens with an intraocular
54 lens (IOL) remains the main goal of cataract surgery. However, modern-day
55 patients are in general more demanding and have higher expectations
56 regarding visual quality, comfort and spectacle independence after IOL
57 implantation. For that reason, cataract surgery is nowadays performed at
58 increasingly earlier age and has become a consolidated option within the
59 portfolio of refractive procedures.¹ This has prompted the constant evolution of
60 IOLs' designs intending to achieve the best visual function possible, especially
61 at intermediate and near distances, while maintaining perceived good image
62 quality at a far distance. Designs currently available on the market are
63 diffractive, refractive or combined refractive/diffractive ones, with low,
64 intermediate or high addition power each providing distinct vision at different
65 distances. With regard to the IOLs' foci feature, the lenses are commonly
66 categorized as multifocal IOLs (i.e., bifocal and trifocal) and extended range of
67 vision, the latter being commonly referred as extended depth of focus (EDOF)
68 IOLs.^{2,3} Effective extension of the depth of focus from distance to intermediate
69 and near distances has been achieved with either, diffractive-based bifocal IOLs
70 which combine low addition and chromatic aberration correction,^{4,5,6} or more
71 recently by means of a refractive-based IOL with alternate zones of different
72 focus power and spherical aberration (SA).⁷ Examples of these EDOF IOLs are
73 Tecnis® Symfony (Johnson & Johnson Vision, Inc., Ireland) and Mini WELL
74 (SIFI, Catania, Italy) respectively.

75 The Tecnis® Eyhance IOL, model ICB00 (ICB-IOL) (Johnson & Johnson
76 Surgical Vision, Inc.) is a new monofocal refractive lens aimed at extending

77 depth of focus in comparison to a standard monofocal IOL. The manufacturer's
78 goal with this new design is to offer the patient better visual acuity at
79 intermediate viewing distances -that is required for many important daily tasks-
80 while maintaining the quality and amount of vision the patient gets for far vision.
81 In addition, and according to the manufacturer, the ICB-IOL should not produce
82 more photic nuisance than a conventional monofocal IOL does,⁸ although this
83 feature has yet to be confirmed by clinical studies. The ICB-IOL incorporates a
84 modified aspheric anterior surface that differs from that of its predecessor, the
85 Tecnis® 1-piece, model ZCB00 (ZCB-IOL), an IOL that has been widely
86 implanted throughout the world and yields well-known outcomes.⁹⁻¹²
87 The aim of this paper is to evaluate in-vitro the optical performance of the new
88 ICB-IOL by comparison to the standard monofocal ZCB-IOL.

89 **METHODS**

90 **Intraocular lenses**

91 Two monofocal IOLs produced by the same manufacturer (Johnson & Johnson
92 Vision, Inc) were included in this study: Tecnis®-1 model ZCB00 (ZCB-IOL) and
93 Tecnis® Eyhance, model ICB00 (ICB-IOL). Both lenses share the same
94 platform, have a biconvex design and are made of the same ultraviolet-light
95 absorbing hydrophobic material with a refractive index of 1.47 (at 35°). In
96 addition, for a 6mm eye entrance pupil (5.3mm at the IOL plane),¹³ they both
97 produce negative 4th-order spherical aberration of -0.27µm.^{8,14} The studied
98 lenses had the same refractive power (20 D).

99 The ZCB-IOL is a standard monofocal lens with an anterior aspheric surface
100 and a posterior spherical one. The optical and clinical performance of this lens
101 have been extensively reported in previous works.^{9,10,12,14-17}.

102 The new ICB-IOL features a modified higher-order aspheric anterior surface
103 intended to produce a continuous power increase from the periphery to the
104 center of the lens. More concretely and according to a manufacture's specialist,
105 whereas power in the ZCB-IOL increases from the periphery to the center of the
106 lens, the power change in the ICB-IOL is continuous, but faster, with most of the
107 change occurring in the central part of the lens.⁸ The posterior surface of the
108 lens is spherical.

109 **Optical quality and halo assessment**

110 The optical performance of the IOLs was evaluated with a test bench that has
111 been described in detail elsewhere,^{17,18} and mainly consists of three parts: the
112 illumination system, the model eye and the image acquisition system (Figure A).
113 Since the ICB and ZCB IOLs share the same hydrophobic acrylic material of the
114 Tecnis® 1-piece family of IOLs, their chromatic properties and spectral
115 performances should be very similar, and then, we have considered green
116 illumination (530nm±20nm) exclusively in our experimental tests. The green
117 LED source illuminated either a four-slit test or a pinhole object for MTF
118 measurements^{17,19} and halo assessment^{20,21,22} respectively. The model eye
119 was formed by an artificial cornea and a wet cell with balanced salt solution
120 where the IOLs were placed. A variable aperture diaphragm, placed in front of
121 the artificial cornea, was used as the entrance pupil (EP) to control the size of
122 the beam on the artificial cornea and hence, the level of corneal SA of the

123 wavefront that impinged upon the tested IOL (Fig. A). Additionally, the EP size
124 also determined the beam size on the IOL plane (referred hereafter to as IOL-
125 pupil).¹³ The ratio IOL-pupil to EP was experimentally calibrated to be 0.56.
126 From now on, all the pupil diameters are referred to the IOL plane.^{18,23} The
127 cornea was an achromatic doublet (Lambda-X, Belgium) that induced +0.175
128 µm of 4th-order SA for a 5.0mm IOL-pupil. The model eye with the IOL formed
129 an image of the test object at its best focus that was projected through a 10X
130 infinity corrected microscope onto an 8-bit CCD camera. All the optical elements
131 in the setup were mounted in high-precision mechanical holders with three axis
132 (X, Y and Z) micrometer-precision adjustments.

133 The modulation transfer function (MTF) of the IOLs placed in the model eye was
134 measured at their best focus plane for distance vision. This focus plane was
135 experimentally determined as the one that maximized the MTF for a 3.0mm
136 IOL-pupil and was set as the origin for defocus (i.e., 0.00D).

137 The through focus MTF curves were obtained between -3.00D to +1.00D in
138 0.10D steps with three IOL-pupil sizes: 2.0mm, 3.0mm and 4.5mm, the last two
139 simulating photopic and mesopic illumination conditions in the clinic.
140 Additionally, the optical quality was also evaluated with the area under the MTF
141 metric (MTFa) given its potential significance as preclinical metric.^{14,17} The
142 MTFa was obtained by integrating the corresponding MTF values from 0 to 50
143 cycles/mm as reported elsewhere.¹⁴ The MTF was computed from the images
144 of the four slit object, and more specifically, from the modulus of the Fourier
145 transform of the line spread function of each slit (i.e., four MTF curves).¹⁹ The
146 mean and standard deviation of the MTF and MTFa were derived from this four

147 measurements. The higher the MTF_a value, the better the optical quality of the
148 IOL.

149 For the halo assessment, we determined the halo energy as illustrated in Figure
150 B. The image provided by the CCD camera (linear scale of intensity), consisted
151 of the sharp and intense image of the pinhole (referred from now on as core)
152 surrounded by a faint halo (Figure BA). When the image was displayed in
153 logarithmic scale of intensity (Fig. BB),^{16, 21, 22, 24} which is a closer representation
154 of how the human eye would see the image, the halo became quite evident.

155 In the image in logarithmic scale, we computed the log-transform energy of the
156 region of interest as:^{16, 18, 25}

$$E_R = \sum_{n \in R} \log(e(n)) \quad \text{Eq.}$$

158 (1)

159 where R stands for either the total image, or the core region - inside the dash
160 black circle in Fig. BB-, or the halo –outside the dash black circle in Fig. BB- ($R=$
161 total, core or halo), n is a pixel contained in the R region, and $e(n)$ is the pixel
162 gray level. For each pupil size, the log-transform energy obtained with Eq. 1 in
163 the core and halo regions (E_{core} and E_{halo} respectively) were compared to E_{total}
164 and expressed as percentages. This normalization is necessary for quantitative
165 comparison of the pinhole images recorded with different pupil sizes and thus
166 with different energy. Moreover, since the human eye responds to differences of
167 energy, we have also computed the non-normalized differences between the
168 core energy and the halo energy, which estimates the weight of the halo in the
169 image: the larger the difference of energy the lower the weight of the halo. The
170 uncertainty in the computed values of the energy was basically due to the
171 precision in the determination of the size of the core. Assuming an uncertainty

172 of ± 1 pixel in the diameter of this region of interest, the highest error
173 corresponded to the lower IOL-pupil (2.0mm) and was 5% for both IOLs.

174 **Wavefront aberrations measurement**

175 The high-order wavefront aberrations of the IOLs were measured from IOL-
176 pupils ranging from 2.0mm to 5.0mm as reported in Reference 26, modifying
177 the optical configuration of the test bench. As shown in the layout of Fig. C, the
178 artificial cornea was removed from the setup and the microscope and CCD
179 camera were replaced by an aberration free collimating lens and a Shack-
180 Hartmann wavefront sensor (HASO 76, Imagine Optics, France). This sensor
181 has an array of 76x100 microlenses, thus providing an excellent spatial
182 resolution during wavefront sampling, and a maximum aperture size of 8.7x11.4
183 mm. Wavefront fitting was made with a linear combination of 32 Zernike
184 polynomials from 3rd to 6th order. Each wavefront was measured three times to
185 obtain the mean value and standard deviation of the Zernike coefficients.

186 Finally, we also computed the wavefront aberrations of the achromatic lens
187 used as artificial cornea. They were obtained versus IOL-pupil size by ray
188 tracing simulation using a dedicated software (Zemax OpticStudio, Zemax
189 Europe Ltd) and the lens parameters (curvature radius, thickness and refraction
190 index) provided by the manufacturer.

191 **RESULTS**

192 **MTF measurements**

193 Figure 1 shows the influence of pupil size on the MTFs of both IOLs. For IOL-
194 pupil sizes of 2.0mm and 3.0mm, the MTF curves of the ZCB-IOL were nearly
195 diffraction limited, while the ones of ICB-IOL were lower, indicating worse

196 optical quality. For larger pupils, the MTF curves of both IOLs tended to
197 decrease and get closer. Overall for both IOLs, it is worth remarking the close
198 coincidence of their curves in the range of spatial frequencies of primary interest
199 (0 to 50 cycles/mm).^{14,17}

200 The MTFa metric versus IOL-pupil size is shown in Figure 2. For pupil sizes
201 lower than 3.5mm, the ICB-IOL had smaller MTFa than ZCB-IOL. For larger
202 pupils, the MTFa values of both IOLs were similar within the experimental
203 uncertainty, and tended to decrease, the later showing the deleterious influence
204 that the increase of the pupil size has on optical quality.

205 The through focus MTFa curves of the ZCB-IOL and ICB-IOL, obtained with
206 IOL-pupils of 2.0mm, 3.0mm and 4.5mm, are shown in Figure 3. Given the
207 monofocal design of both lenses, the curves showed just one peak of maximum
208 MTFa that corresponds to the best focus of the lenses for distance vision.

209 For both IOLs, the smaller the pupil, the wider the MTFa peak, proving that
210 regardless of the IOL design, there was an effect of focus extension produced
211 as a consequence of reducing the pupil size. With a 2.0mm IOL-pupil, the
212 maximum MTFa of the ZCB-IOL was higher than the one of the ICB-IOL
213 (45.69±0.23 versus 40.61±0.49) and very interestingly, there was a myopic shift
214 of -0.40D in the position of the MTFa peak of the ICB-IOL (Fig. 3A). Even at this
215 position, however, the MTFa value reached by the ICB-IOL was not higher than
216 that of ZCB-IOL. With a 3.0 mm IOL-pupil (Fig. 3B), the maximum MTFa of the
217 ZCB-IOL was still slightly higher (46.11±0.58 versus 42.79±0.77 in the case of
218 the ICB_IOL) while such differences between both IOLs practically vanished

219 with the 4.5mm IOL-pupil (Fig. 3C). With this pupil, the MTFa peak of both IOLs
220 shifted slightly towards hyperopic defocus (+0.20D).

221 **Wavefront Aberration**

222 The cornea of our model eye induced positive 4th-order spherical aberration
223 (SA) (Figure 4) to mimic the natural aberration of the human cornea. The
224 maximum SA was +0.175 μm for a 5.0mm pupil at the IOL plane. Although this
225 value is somehow lesser than the amount reported on average for the human
226 cornea, +0.27 μm for a 6.0mm entrance pupil (5.3mm at the IOL plane),²⁷ Wang
227 et al.,²⁸ found that 15.4% of their patients had corneas with SA values smaller
228 than +0.2 μm .

229 The most significant high-order aberration found with both IOLs versus pupil
230 size was negative 4th-order SA (Figure 4), logically intended to compensate for
231 the positive corneal SA. With the 3.0mm IOL-pupil, the wavefront aberration of
232 the ICB-IOL also showed small contributions of positive 6th-order (0.028 ± 0.001)
233 and negative 8th-order SA (-0.018 ± 0.001). The rest of high-order aberration
234 terms were negligible small for all pupil sizes.

235 For IOL-pupil sizes lower than 3.5mm, the ICB-IOL lens had SA values more
236 negatives than ZCB-IOL. For instance, with a 2.0mm IOL-pupil the SA of ICB-
237 IOL is, in absolute value, 3.7 times larger than the SA of ZCB-IOL (-
238 $0.056 \pm 0.003 \mu\text{m}$ versus $-0.015 \pm 0.003 \mu\text{m}$ respectively). For larger pupils both
239 lenses had very similar SA values.

240 **Halo assessment**

241 The images of the pinhole as a function of the pupil size, obtained with the
242 model eye including either ZCB-IOL or ICB-IOL, are shown in logarithmic scale
243 in Figure 5 since it better approaches human perception. With both lenses,
244 halos around the pinhole image could be hardly observed for IOL-pupil sizes up
245 to 3.0mm, but became apparent for larger IOL-pupil sizes of 4.0mm and 5.0mm.

246 The relative energy of the core and halo to total energy (all calculated with Eq.
247 1), are shown versus IOL-pupil size in Figures 6(A) and 6(B) respectively. On
248 the other hand, Figure 6(C) shows the non-normalized energy difference
249 between the core and halo regions.

250 In the case of the ZCB-IOL and for IOL-pupils ranging from 2.0mm up to 3.5
251 mm, a constant, high fraction of the energy ($\approx 60\%$) was correctly focused on the
252 core, while with the ICB-IOL we found less energy correctly focused.
253 Conversely, there was more energy spread to the halo with the ICB-IOL. Not
254 surprisingly, ICB-IOL showed in this pupil range, smaller values of the energy
255 difference than ZCB-IOL (Fig. 6(C)), indicating images with more significant
256 halos with ICB-IOL. This trend could already be acknowledged from Fig. 5
257 where the halo for the 2.0mm IOL-pupil with the ICB-IOL was clearly more
258 visible and larger than that of ZCB-IOL, but was only slightly more visible and
259 larger for the 3.0 mm IOL-pupil. The differences between the two IOLs tended
260 to reduce for increasing pupils. As such, for mesopic and scotopic pupils
261 ($\geq 4.0\text{mm}$), both IOLs exhibited very close results. The halo energy for both IOLs
262 increased with pupil size (Fig. 6(B)), reaching similar maximum values for the
263 5.0mm IOL-pupil: $71.9 \pm 1.8\%$ (ZCB-IOL) and $68.5 \pm 1.7\%$ (ICB-IOL). Closely
264 related, the smaller value of energy difference occurred for this IOL-pupil (Fig.
265 6(C)), which accounts for the significant halos observed in Fig. 5 with both IOLs.

266 A summary of the results versus IOL-pupil size is presented in Table 1.

267 **DISCUSSION**

268 To our knowledge, this is the first paper that studied *in-vitro* the optical
269 performance of the new monofocal Tecnis® Eyhance model ICB00 (ICB-IOL)
270 whose optical design aims at extend the DOF in comparison to a standard
271 monofocal lens, to provide better intermediate vision, while keeping similar
272 distance vision and comparable incidence of photic phenomena (glare and
273 halo). A meaningful comparison has been carried out by choosing the Tecnis®
274 1-piece model ZCB00 (ZCB-IOL) as the standard monofocal lens to compare
275 with, because both IOLs, manufactured by the same company, share basic
276 features such as platform and material. Therefore, the fundamental difference
277 between both designs is the modified aspheric anterior surface of the ICB-IOL
278 that is referred by the manufacturer as a continuous higher-order aspheric
279 surface. We have found that this modification in the optical design of the ICB-
280 IOL has a measurable impact on the optical quality (Figures 1, 2 and 3), the SA
281 (Fig. 4) and halo energy (Fig. 6) for relatively small IOL-pupil sizes (below
282 3.5mm). In contrast, for larger IOL-pupils (equal and above of 3.5mm) the
283 results of the new ICB-IOL tend to be very similar to the standard monofocal
284 ZCB-IOL.

285 The negative values of SA versus pupil obtained with the ZCB-IOL (Fig. 4) are
286 in good agreement with previous results reported with multifocal and EDOF
287 Tecnis® IOLs that share the same aspheric design.^{26,29,30} The measured
288 experimental values of $-0.28 \pm 0.01 \mu\text{m}$ (ZCB-IOL) and $-0.27 \pm 0.01 \mu\text{m}$ (ICB-IOL)
289 for 5.0mm IOL-pupil, would fully compensate for the $+0.27 \mu\text{m}$ value of the

290 corneal SA of a representative average human cornea.^{27,28} Interestingly, we
291 have found differences between the SA values of both IOLs for IOL-pupils lower
292 than 3.5mm and it is in this range of small pupils where worse MTF curves (Fig.
293 1) and lower MTF_a values (Fig. 2) are obtained for the new ICB-IOL in
294 comparison to ZCB-IOL. Moreover, comparing the through focus MTF_a curves
295 of the two IOLs (Fig. 3), the largest differences occurred with the smallest IOL-
296 pupil (2.0mm). For this pupil in particular, the MTF_a curves of both IOLs are
297 considerably broader, proving that small pupils are an effective strategy to
298 expand DOF in general,³¹ although they require good lighting. More importantly,
299 in comparison to the standard ZCB-IOL, the MTF_a curve of the new ICB-IOL
300 showed a myopic shift of -0.40D (Fig.3A). To explain this result, we recall that
301 for small pupils the ICB-IOL induced more negative SA than the standard ZCB-
302 IOL, and they were larger (in absolute value) than the SA of the cornea (Figure
303 4). Then, the model eye with the ICB-IOL must have a remaining negative SA
304 as a result of the insufficient compensation between the positive and the
305 negative SA values of the cornea and ICB-IOL respectively. In a converging
306 optical system with negative SA, the paraxial rays have more dioptric power
307 than the peripheral rays. Thus, in eyes with relatively large pupil (e.g. mesopic
308 illumination conditions) and negative SA, the emmetropia condition is achieved
309 when the circle of least confusion lies on the retina, with the paraxial and
310 peripheral rays focused in front of and behind the retina respectively. Since the
311 eye focusing with small pupils relays only in the paraxial rays, there would be a
312 myopic shift of the best focus condition, as experimentally observed in the case
313 of the ICB-IOL (Fig. 3A). This result confirms the power increase from the
314 periphery to the center in the design of the new ICB-IOL. More concretely and in

315 comparison to the standard ZCB-IOL, there is an additional power close to
316 +0.50D in the central 2.0mm region of the lens, which is based on the larger
317 negative SA of the ICB-IOL with this IOL-pupil (Figure 4). Since the differences
318 between the SA of the two IOLs decreased from 2.0mm to 3.0mm and were
319 practically equal for IOL-pupils larger than 3.5mm, one can conclude that the
320 aspheric curvature of the ICB-IOL originates the +0.50D additional power in the
321 central 2.0 to 3.0mm region and decreases towards the periphery of the lens.
322 This is the basis for the intermediate performance and extension of the depth of
323 focus with this new design, which could improve intermediate vision, especially
324 in high light conditions and/or small pupils.

325 Several studies have additionally shown that higher-order aberrations,
326 particularly SA, helps to increase the DOF.^{32,33} However, the addition of SA to
327 increase the DOF has the potential drawback of lowering the visual acuity
328 (VA).³⁴ In the clinic, Rocha et al.,³⁵ and Marcos et al.,³⁴ found larger DOF in
329 patients implanted with spherical IOLs (i.e., eyes with higher SA) than in
330 patients with aspheric IOLs that reduced the total SA of the eye. Other studies
331 however, failed to find statistically significant differences in DOF between
332 patients implanted with aspheric IOLs with negative SA, aspheric aberration-
333 free IOLs, and spherical IOLs.^{36,37} Neither did XianHui et al. find differences in
334 the DOF of eyes with different amounts of 4th-order corneal SA implanted with
335 the same aberration free IOL model.³⁸ More recently, Belluci et al.³⁹ reported, in
336 comparison to an aspheric monofocal IOL, greater depth of focus with the Mini
337 Well® IOL, a lens based on alternating positive and negative SA in the central
338 3.0mm optical zone. Camps et al.²⁶ measured large negative 4th-order (-
339 0.13±0.01 μm) and positive 6th-order (0.12±0.01 μm) SA within this 3.0mm

340 zone of the Mini Well IOL. Since the experimental SA values of the ICB-IOL are
341 much smaller (Fig. 4), it can be ruled out that the DOF expansion with this new
342 monofocal lens be based on a SA design.

343 With regard to the differences on optical quality between the ZCB-IOL and ICB-
344 IOL accounted by the MTFa metric and its implication in the clinic, Alarcon et
345 al.¹⁴ found a high correlation between the MTFa metric and clinical VA of
346 pseudophakic patients. The results led the authors to suggest that the MTFa
347 metric could predict clinical average VA, thus becoming a preclinical metrics.
348 More recently, Vega et al.¹⁷ stated that the estimation of achievable VA, as non-
349 linear function of variable MTFa, showed limiting behavior for IOLs with larger
350 MTFa values, i.e. lenses with higher imaging quality. As a consequence,
351 beyond an MTFa threshold, VA tended asymptotically to the best value clinically
352 achieved in the patients, and any further increase in the imaging quality of an
353 IOL (i.e., MTFa values above the threshold) did not translate into VA
354 improvement. With the 3.0mm IOL-pupil, the best MTFa of both IOLs (ICB-
355 IOL=43.15±0.43, ZCB-IOL=46.10±0.33) were far larger than the reported
356 threshold for the MTFa (≈ 20),¹⁷ and then, one would not expect that the best VA
357 with the new ICB-IOL model be worse than with the standard ZCB-IOL.
358 Nevertheless, clinical studies are needed to either confirm or refute this
359 prediction.

360 Concerning the halos, the ICB-IOL lens showed on the optical bench, for
361 mesopic and scotopic pupils ($\geq 4.0\text{mm}$), similar values of energy correctly
362 focused on the core region (or conversely spoiled in the halo) compared to the
363 standard monofocal (Fig. 6). Thus, the new design of the ICB-IOL is more likely
364 to induce similar level of photic phenomena, if any, as to the ZCB-IOL, although

365 clinical studies, yet to be performed, are mandatory to confirm or refute this
366 prediction.

367 Regarding the differences found in the spherical aberration (SA) within the small
368 central area between the two IOLs, Taketani and Hara,⁴⁰ have shown that SA
369 was negatively correlated with dioptric power in the case of the Tecnis ZA9003
370 aspheric IOL (a lens made of the same material and with the same optical
371 design as the Tecnis ZCB00 IOL of our study). Then, additional work is still
372 necessary to check if the differences of SA we have found between ICB00 and
373 ZCB00 with lenses of +20.0D, are maintained for other dioptric powers. This
374 would be especially relevant in the case of low and high dioptric values because
375 it could provide very valuable information on the performance of the new ICB00
376 design in the case of highly myopic and hyperopic eyes respectively.

377 Finally, our results have shown that the ICB-IOL is a modified monofocal lens
378 with 0.50D of additional power in the central 2 mm zone. With regard to the
379 potential impact that this new design may have on near vision, one could argue
380 that pupil miosis would occur when looking at near, a situation under which the
381 ICB-IOL has shown the capability, in the optical bench, of producing a myopic
382 shift related to the power increase on the central region of the lens. However,
383 near tasks at 30-40 cm that demand good quality of vision require add powers
384 of 3.0-2.5D, which are significantly larger than the maximum additional power
385 (\approx 0.50D) that the lens is able to provide. Furthermore, despite the potential to
386 extend the range of vision of the ICB-IOL, it is yet a monofocal design.
387 Preliminary clinical results (N. Garzón, private communication) evidences that
388 the ICB-IOL defocus curve is slightly broader than the one of standard ZCB-
389 IOL, but still has the typical ‘Λ-like monofocal shape’ with a single visual acuity

390 (VA) peak for distance vision (0.0D defocus). Thus, it would not be realistic to
391 expect that the ICB-IOL could compete in near vision with multifocal IOLs that
392 provide the pseudophakic patient with additional VA peak at near

393 To summarize, according to our findings the design strategy of the new ICB-IOL
394 to extend DOF is based on a continuous power increase towards the center of
395 the lens as a result of the increased amount of negative SA that occurs in the
396 central region of the lens (IOL-pupils lower than 3.5mm). The halos measured
397 on the optical bench are comparable to a standard monofocal IOL.

398 **References**

- 399 1. Wolffsohn JS, Davies LN. Presbyopia: Effectiveness of correction
400 strategies. *Prog Retin Eye Res.* 2019;68:124-143.
- 401 2. Breyer DRH, Kaymak H, Ax T, Kretz FTA, Auffarth GU, Hagen PR.
402 Multifocal Intraocular Lenses and Extended Depth of Focus Intraocular
403 Lenses. *Asia-Pacific J Ophthalmol.* 2017;6(4):339-349.
- 404 3. Kondylis G, Klavdianou O, Palioura S. Multifocal and extended depth of
405 focus intraocular lenses. *Ann Eye Sci.* 2019;4:5.
406 doi:10.21037/aes.2019.01.01
- 407 4. Weeber HA, Meijer ST, Piers PA. Extending the range of vision using
408 diffractive intraocular lens technology. *J Cataract Refract Surg.*
409 2015;41(12):2746-2754.
- 410 5. Millán MS, Vega F. Extended depth of focus intraocular lens : chromatic
411 performance. *Biomed Opt Express.* 2017;8(9):4294-4309.
- 412 6. Łabuz G, Papadatou E, Khoramnia R, Auffarth GU. Longitudinal
413 Chromatic Aberration and Polychromatic Image Quality Metrics of
414 Intraocular Lenses. *J Refract Surg.* 2018;34(12):832-838.
- 415 7. Bellucci R, Curatolo MC. A New Extended Depth of Focus Intraocular
416 Lens Based on Spherical Aberration. *J Refract Surg.* 2017;33(6):389-394.
- 417 8. Alarcon A, Koopman B, Canovas C, Domingo J, Piers P. Optical and
418 predicted visual performance of the next generation Tecnis monofocal
419 intraocular lens. Presented at the 23rd ESCRS winter meeting; Athens.
420 February 15-17; 2019.

- 421 9. Yamauchi T, Tabuchi H, Takase K, Ohsugi H, Ohara Z, Kiuchi Y.
- 422 Comparison of Visual Performance of Multifocal Intraocular Lenses with
- 423 Same Material Monofocal Intraocular Lenses. *PLoS One*. 2013;8(6):1-7.
- 424 10. Pedrotti E, Bruni E, Bonacci E, Badalamenti R, Mastropasqua R, Marchini
- 425 G. Comparative Analysis of the Clinical Outcomes With a Monofocal and
- 426 an Extended Range of Vision Intraocular Lens. *J Refract Surg*.
- 427 2016;32(7):436-442.
- 428 11. Zhao Y, Wang Z, Tian X, Wang X, Gao X. Comparative study of visual
- 429 function and ocular aberrations of two different one-piece designed
- 430 hydrophilic acrylic intraocular lens. *Int Ophthalmol*. 2018;38(3):1169-1175.
- 431 12. Altemir-Gomez I, Millán MS, Vega F, et al. Comparison of visual and
- 432 optical quality of monofocal versus multifocal intraocular lenses. *Eur J*
- 433 *Ophthalmol*. 2019. doi:10.1177/1120672119827858
- 434 13. Norrby S, Piers P, Campbell C, van der Mooren M. Model eyes for
- 435 evaluation of intraocular lenses. *Appl Opt*. 2007;46(26):6595-6605.
- 436 14. Alarcon A, Canovas C, Rosen R, et al. Preclinical metrics to predict
- 437 through-focus visual acuity for pseudophakic patients. *Biomed Opt*
- 438 *Express*. 2016;7(5):1877-1888.
- 439 15. Kretz FT, Tandogan T, Khoramnia R AG. High order aberration and
- 440 straylight evaluation after cataract surgery with implantation of an
- 441 aspheric, aberration correcting monofocal intraocular lens. *Int J*
- 442 *Ophthalmol*. 2015;18;8(4):736-741.
- 443 16. Yoo Y-S, Whang W-J, Byun Y-S, et al. Through-Focus Optical Bench

- 444 Performance of Extended Depth-of-Focus and Bifocal Intraocular Lenses
445 Compared to a Monofocal Lens. *J Refract Surg.* 2018;34(4):236-243.
- 446 17. Vega F, Millán MS, Garzón N, Altemir I, Poyales F, Larrosa JM. Visual
447 acuity of pseudophakic patients predicted from in-vitro measurements of
448 intraocular lenses with different design. *Biomed Opt Express.*
449 2018;9(10):4893-4906.
- 450 18. Vega F, Alba-Bueno F, Millán MS. Energy distribution between distance
451 and near images in apodized diffractive multifocal intraocular lenses.
452 *Invest Ophthalmol Vis Sci.* 2011;52(8):5695-5701.
- 453 19. Simpson MJ. Diffractive multifocal intraocular lens image quality. *Appl
454 Opt.* 1992;31(19):3621-3626.
- 455 20. Vega F, Alba-Bueno F, Millán MS, Varon C, Gil MA, Buil JA. Halo and
456 Through-Focus Performance of Four Diffractive Multifocal Intraocular
457 Lenses. *Investig Ophthalmology Vis Sci.* 2015;56(6):3967-3975.
- 458 21. Alba-Bueno F, Vega F, Millan MS. Halos and multifocal intraocular
459 lenses: Origin and interpretation. *Arch Soc Española Oftalmol (English
460 Ed).* 2014;89(10):397-404.
- 461 22. Alba-Bueno F, Garzón N, Vega F, Poyales F, Millán MS. Patient-
462 Perceived and Laboratory-Measured Halos Associated with Diffractive
463 Bifocal and Trifocal Intraocular Lenses. *Curr Eye Res.* 2018;43(1):35-42.
- 464 23. Vega F, Alba-Bueno F, Millán MS. Energy efficiency of a new trifocal
465 intraocular lens. *J Eur Opt Soc Rapid Publ.* 2014;9:14002.
466 doi:10.2971/jeos.2014.14002

- 467 24. Alba-Bueno F, Millán MS, Vega F. Optical Characterization of Intraocular
468 Lenses. *Opt. Pura Apl.* 2017;50(1): 63-73
- 469 25. Vega F, Millán MS, Vila-Terricabras N, Alba-Bueno F. Visible Versus
470 Near-Infrared Optical Performance of Diffractive Multifocal Intraocular
471 Lenses. *Investig Ophthalmology Vis Sci.* 2015;56(12):7345-7351.
- 472 26. Camps VJ, Tolosa A, Piñero DP, De Fez D, Caballero MT, Miret JJ. In
473 Vitro Aberrometric Assessment of a Multifocal Intraocular Lens and Two
474 Extended Depth of Focus IOLs. *J Ophthalmol.* 2017; 7 pages.
475 doi:10.1155/2017/7095734
- 476 27. Guirao A, Redondo M, Artal P. Optical aberrations of the human cornea
477 as a function of age. *J Opt Soc Am A Opt Image Sci Vis.*
478 2000;17(10):1697-1702.
- 479 28. Wang L, Dai E, Koch DD, Nathoo A. Optical aberrations of the human
480 anterior cornea. *J Cataract Refract Surg.* 2003;29(8):1514-1521.
- 481 29. Kim MJ, Zheleznyak L, Macrae S, Tchah H, Yoon G. Objective evaluation
482 of through-focus optical performance of presbyopia-correcting intraocular
483 lenses using an optical bench system. *J Cataract Refract Surg.*
484 2011;37(7):1305-1312.
- 485 30. Gatinel D, Loicq J. Clinically Relevant Optical Properties of Bifocal,
486 Trifocal, and Extended Depth of Focus Intraocular Lenses. *J Refract
487 Surg.* 2016;32(4):273-280.
- 488 31. Xu R, Wang H, Jaskulski M, Kollbaum P, Bradley A. Small-pupil versus
489 multifocal strategies for expanding depth of focus of presbyopic eyes. *J*

- 490 *Cataract Refract Surg.* 2019;45(5):647-655.
- 491 32. Rocha KM, Vabre L, Chateau N, Krueger RR. Expanding depth of focus
492 by modifying higher-order aberrations induced by an adaptive optics
493 visual simulator. *J Cataract Refract Surg.* 2009;35(11):1885-1892.
- 494 33. Benard Y, Lopez-Gil N, Legras R. Subjective depth of field in presence of
495 4th-order and 6th-order Zernike spherical aberration using adaptive optics
496 technology. *J Cataract Refract Surg.* 2010;36(12):2129-2138.
- 497 34. Marcos S, Barbero S, Jiménez-Alfaro I. Optical quality and depth-of-field
498 of eyes implanted with spherical and aspheric intraocular lenses. *J*
499 *Refract Surg.* 2005;21(3):223–235.
- 500 35. Rocha KM, Soriano ES, Chamon W, Chalita MR, Nosé W. Spherical
501 Aberration and Depth of Focus in Eyes Implanted with Aspheric and
502 Spherical Intraocular Lenses. A Prospective Randomized Study.
503 *Ophthalmology.* 2007;114(11):2050-2054.
- 504 36. Santhiago MR, Netto M V., Barreto JR. J, et al. Wavefront analysis,
505 contrast sensitivity, and depth of focus after cataract surgery with
506 aspherical intraocular lens implantation. *Am J Ophthalmol.*
507 2010;149(3):383-389.
- 508 37. Shentu X, Tang X, Yao K. Spherical aberration, visual performance and
509 pseudoaccommodation of eyes implanted with different aspheric
510 intraocular lens. *Clin Experiment Ophthalmol.* 2008;36(7):620-624.
- 511 38. Gong X-H, Zheng Q-X, Wang N, et al. Visual and optical performance of
512 eyes with different corneal spherical aberration implanted with aspheric

- 513 intraocular lens. *Int J Ophthalmol.* 2012;5(3):323-328.
- 514 39. Bellucci R, Cargnoni M, Bellucci C. Clinical and aberrometric evaluation of
515 a new extended depth-of-focus intraocular lens based on spherical
516 aberration. *J Cataract Refract Surg.* 2019: 8 pages. In Press, Corrected
517 Proof, Available online 24 May 2019 doi:10.1016/j.jcrs.2019.02.023.
- 518 40. Taketani F, Hara Y. Characteristics of spherical aberrations in 3 aspheric
519 intraocular lens models measured in a model eye. *J Cataract Refract
520 Surg.* 2011;37(5):931–936.

521

522 **Figure Captions**

523 Figure 1. MTF curves for IOL-pupil sizes ranging from 2.0mm to 5.0mm at the
524 best focus of the (—) ZCB-IOL and (—) ICB-IOL. The inserts show the MTF
525 curves in the range of spatial frequencies (0 to 50 cycles/mm) used to compute
526 the MTF_a metric.

527 Figure 2. MTF_a values (average± standard deviation) for IOL-pupil sizes ranging
528 from 2.0mm up to 5.0mm of the (●) ZCB-IOL and (●) ICB-IOL. These values
529 were obtained upon integration of the corresponding MTF curves between 0 to
530 50 cycles/mm.

531 Figure 3. Through focus MTF_a curves of the (—) ZCB-IOL and (—) ICB-IOL
532 obtained with IOL-pupil sizes of (A) 2.0mm, (B) 3.0mm and (C) 4.5mm. The
533 arrows indicate the myopic (A) and hyperopic (B) shifts of the MTF_a peaks.

534 Figure 4. 4th-order SA versus IOL-pupil size obtained separately for the artificial
535 cornea (black bars), and each IOL: ZCB-IOL (red bars) and ICB-IOL (blue bars).

536 Figure 5. Images of the pinhole object formed by the model eye including either
537 ZCB-IOL or ICB-IOL at their best focus for increasing pupil sizes. The images
538 are displayed in logarithmic scale of intensity.

539 Figure 6. Relative core (A) and halo energy (B) to total energy (all calculated
540 with Eq. 1) and (C) non-normalized energy difference between core and halo,
541 versus IOL-pupil size, obtained with the ZCB-IOL (red bars) and ICB-IOL (blue
542 bars).